Alphaplast

DMG Chemisch-Pharmazeutische Fabrik GmbH
Elbgastrasse 248 - 22547 Hamburg - www.dmg-dental.com
091833/#0.1_2013-09

0482

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Deutsch

Produktbeschreibung

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English

Product description

Traditional self-curing quartz composite.

Indications

Class III, IV and V restorations.

Contra-indications

- Applying Alphaplast is contra-indicated if a dry working area or the specified application technique is not possible.
- Do not use directly on exposed pulp or in cases of allergies to any of the ingredients.

Basic Safety Instructions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of inadvertent skin contact, wash the affected area thoroughly with soap and water and seek medical attention, if necessary.
- Avoid contact with the eyes! In the event of inadvertent contact with eyes, immediately rinse the eyes thoroughly with plenty of water and consult a physician if necessary.

Side effects

There are no known side effects to date. Isolated cases of hypersensitivity to components of the material cannot be excluded.

Interactions with other materials

- Products containing eugenol, moisture and air containing oil can affect the curing of Alphaplast. Avoid products containing eugenol, moisture and air containing oil.
- Do not use metal spatulas.
- The use of cationic mouthwashes, plaque indicators and chlorhexidine may lead to discoloration.

Timing

Mixing time	min. 20 sec
Working time (including mixing time)*	1:30 min
Curing time	2:00 to 3:00 min

*The time indicated applies to a room temperature of 23 °C. Higher temperatures shorten this time and lower temperatures increase this time.

Instructions for use

- Wipe excess composite off the filling instrument before it hardens.

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Dosage and mixing

- Mix the paste in a 1:1 ratio. To do this, place the desired amount of base paste with one end of a plastic spatula onto the mixing block and using the other end of the plastic spatula, place the same amount of catalyst paste on the mixing block (in each case approximately half of the volume of the restoration).
- Mix the pastes for at least 20 s until the mixture is uniform.

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