# Honigum Pro

Heavy Scan

## Inhalt | Contents

Gebrauchsinformation – Deutsch	4
Instructions for use – English	13
Mode d'emploi – Français	22
Istruzioni per l'uso – Italiano	31
Instrucciones de empleo – Español	40
Instruções de utilização – Português	49
Gebruiksaanwijzing – Nederlands	58
Brugsanvisning – Dansk	67
Bruksanvisning – Svenska	76
Instrukcja użycia – Polski	85
Инструкция по применению – Русский	94
使用说明书 – 中文版	104
取扱説明書 – 日本語	111
사용 설명서 – 한국어	120

## Produktbeschreibung

---

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---

Vinylpolysiloxane, Siliciumdioxid, Pigmente, Additive, Platinkatalysator

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## Product description

---

Honigum Pro-Heavy Scan is a scannable precision pre-impression material based on addition-cured silicones. The material is suitable for all scanners designed for scanning impression materials (e.g. the D700 dental scanner from 3Shape). Additional surface treatment is not required.

The material is available in the Automix cartridge and in the ready-to-use pre-filled Mix-Star cartridge, which is suitable for use in the DMG MixStar-eMotion unit and other commercially available mixing machines.

## Indications

---

- Crown and bridge impressions
- Inlay and onlay impressions
- All types of pickup impressions, e.g. for implants and telescopes

## Basic safety information

---

- For dental use only!
- Keep out of the reach of children!

## Instructions for use

---


- The specified working times (see "Technical data") apply at a room temperature of 23 °C/72 °F and a normal relative humidity of 50 %. Higher temperatures shorten these times while lower temperatures increase them.
- Rinse off residual retraction solution thoroughly with water, otherwise this may increase the setting time in some instances.
- Do not allow the impression to come into contact with organic solvents or liquids containing solvents as this may otherwise result in swelling of the material, leading to inaccuracies in the cast.
- Do not allow the impression material to come into contact with latex gloves. DMG

recommends the use of PVC or PE-based gloves.


- Do not combine the material with condensation-curing silicone or polyether impression materials.

### Using the MixStar cartridge in the MixStar-eMotion unit from DMG


- ▶ Follow the instructions for use for the "MixStar-eMotion" unit!


1. Slide the locking mechanism upwards in order to open the cartridge. Do not remove the locking mechanism from the cartridge.

2. Remove the seal from the cartridge and discard. Attach a new mixing tip.

3. Secure the mixing tip by closing the locking mechanism.

Insert the cartridge into the MixStar-eMotion unit. A brief acoustic signal confirms correct insertion of the cartridge.

**Note:** Prior to using for the first time, prepare fresh cartridges by discarding approx. the first 3 cm of material.

After application, leave the used mixing tip on the cartridge as a seal.

## Using the MixStar cartridge in mixing units from other manufacturers

---

- ▶ Follow the instructions for use for the mixing unit used!
- ▶ First insert the cartridge into the mixing unit without the mixing tip. Otherwise, correct mixing of the material cannot be guaranteed!
- ▶ In this case, proceed as described under "Changing the mixing tip"

## Changing the mixing tip

---


- ▶ Follow the instructions for use for the mixing unit used!

**Note:** Attach a new mixing tip each time before mixing.

1. Slide the locking mechanism upwards.
2. Remove the old mixing tip as well as any hardened material blocking the cartridge openings.

- Place the mixing tip securely on the respective cartridge openings.
- Lock the mixing tip in place with the locking mechanism.

## Using the Automix cartridge

### Inserting the cartridge

- Push up the lever [A] at the back of the dispensing gun and pull back the slide [B] completely.
- Lift up the plastic lever [C], insert the cartridge and lock it using the plastic lever [C].

**Note:** Make sure that the notches on the cartridge line up with those on the dispensing gun.

### Attaching a mixing tip

- Remove the cartridge cap [D] or used mixing tip after turning it 90° counterclockwise, and discard it.
- Attach a new mixing tip [E].

**Note:** Make sure that the notches on the mixing tip line up with those on the cartridge.

- Turn the mixing tip clockwise 90° to lock it.

### Applying material

- ▶ When using a cartridge for the first time, squeeze out approx. a pea-sized amount of material and discard it.
- ▶ Mix the material in the mixing tip using the lever [F] on the dispensing gun and apply directly.

**Note:** After application, leave the used mixing tip on the cartridge as a seal!

## Removing the cartridge

1. Push up the lever [A] at the back of the dispensing gun and pull back the slide [B] completely.
2. Lift up the plastic lever [C] and remove the cartridge.

## Recommended use

---

- Honigum Pro-Heavy Scan can be used in all conventional impression trays.
- In order to improve the adhesion of the material to metal or plastic impression trays, a tray adhesive for A-silicones (e.g. tray adhesive, DMG) should be used.
- The final impression can be optically registered using an appropriate scanner or cast in conventional fashion using plaster.

## Double mixing impression

- ▶ The individual times for the materials used are provided in the "Technical data" section. Follow the corresponding instructions for use!

1. Load the impression tray with Honigum Pro-Heavy Scan. When doing so, hold the impression tray at a slight angle beneath the mixing tip.

**Note:** To avoid air bubbles, the end of the mixing tip should remain immersed in the material throughout application.

2. Inject around the preparation with Honigum Pro-Light Scan.

**Note:** Filling of the impression tray and injection should be concluded at the same time, as the higher temperature in the mouth causes the injected material to set more quickly than the material in the impression tray.

3. Position the filled impression tray slowly and gently in the mouth within the working time and form the preparation.
4. Remove the impression from the mouth once the setting time has ended, see „Removal“.

## Final impression

- ▶ The individual times for the materials used are provided in the "Technical data" section. Follow the corresponding instructions for use!

**Note:** Do not allow more than 30 minutes to pass between completion of the pre-impression and the final impression and do not use alcohol for cleaning during this time, as it is otherwise not possible to guarantee a secure bond between the two impression materials.

1. Cut generous relief channels into the pre-impression and remove any imperfections.
2. After testing repositioning, rinse the impression thoroughly with water and dry.
3. Apply Honigum Pro-Light Scan around the preparation and on the pre-impression.
4. Position the filled impression tray slowly and gently in the mouth within the working time, initially exerting brief pressure.
5. Remove the impression from the mouth once the setting time has ended, see „Removal“.

## Removal

---

1. In the area where the teeth have been prepared, remove the impression parallel to the longitudinal axis of the teeth.
2. After removal, cut the impression free as far as possible in order to facilitate optimum scanner access to the relevant areas."
3. Then proceed as described under „Cleaning and disinfection of the final impression“.

## Cleaning and disinfection of the final impression

---

- ▶ Only disinfect the final impression using dipping solutions. DMG recommends the use of a 2% glutaraldehyde solution.
- ▶ Ensure that the disinfectant used is indicated for use with addition-cured silicone impression material.
- ▶ The manufacturer's instructions for the disinfectant solution must be followed exactly!

- ▶ Dry the final impression and store until casting, see "Storing the impression".
- ▶ For more information, see the statement issued by the German Society for Dental and Oral Medicine (DGZMK). Statement by the DGZMK 3/93 V 1.0; [www.dgzmk.de](http://www.dgzmk.de)
- ▶ ADA recommends following Council on Dental Materials, Instruments and Equipment. Disinfection of impressions. JADA 1991;122(3):110; [www.ada.org](http://www.ada.org)

## **Optical scanning**

---

The quality of the optical scanning results may vary depending on the impression scanner. For optimum results, DMG recommends the use of the D700 dental scanner from 3Shape. Observe the manufacturer's instructions for the impression scanner used.

## **Casting and galvanization**

---

- ▶ After removing it from the mouth, wait at least 30 minutes before casting the impression.

DMG recommends using hard plaster (ISO 6873; Type 3, e.g. Moldano, Bayer) for the cast and ultra hard dental plaster for the molds (ISO 6873; Type 4, e.g. Fujirock II, GC). The impression can be copper or silver-plated.

## **Storing the impression**

---

- ▶ Store the impression in a dry place at room temperature (15–25 °C / 59 – 77 °F)!
- ▶ Keep away from direct sunlight!

The impression is dimensionally stable and can be cast without loss of accuracy even after several months.

## Technical data

Honigum Pro-Heavy Scan		
Working time* [min]	≤ 2:45	≤ 2:45
Setting time in the mouth [min]	≥ 3:15	≥ 3:15
Dimensional change after 24 hours [%]	≤ 0.40	≤ 0.40
Residual compression set [%]	≤ 0.50	≤ 0.50
Compression set [%]	≈ 1.9	≈ 1.9

\*The specified times apply at a room temperature of 23 °C and a normal relative humidity of 50 %. Higher temperatures shorten these times while lower temperatures increase them.

The working time describes the period after mixing has begun until the maximum time by which the impression tray must be positioned in the patient's mouth.

## Disinfection of the Automix-Dispenser

The Automix-Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute). The Automix-Dispenser is also autoclavable.

## Composition

Vinyl polysiloxanes, silicon dioxide, pigments, additives, platinum catalyst

## Classification

ISO 4823:2000, Type 1: Heavy-bodied consistency

## Storage and shelf life

---

- ▶ Store in a dry place at room temperature (15–25 °C / 59–77 °F)!
- ▶ Do not store in the fridge or freezer!
- ▶ Keep away from direct sunlight!
- ▶ After application, leave the used mixing tip on the cartridge as a seal!
- ▶ Do not use after the expiration date!

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**DMG** Chemisch-Pharmazeutische Fabrik GmbH