# LuxaCam

## PEEK – Ronde

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## Product description

---

LuxaCam PEEK blanks are tooth-colored, metal-free milling blanks for the fabrication of dental restorations. This material is not radiopaque.

## Indications

---

### Tight-fitting tooth replacement

- Crown and bridge structures in the front and posterior region, 4 elements with up to 2 intermediate elements

### Removable dentures

- Secondary bridge structures
- Telescoping primary and secondary crowns and structures
- One-piece bases on primary sections with cover dentures
- Clamp-retained dentures (model casting)
- Prosthesis bases (e.g. total prostheses)

### Implantology, individual abutments

- Superstructures in the front and posterior region

## Contraindications

---

- Do not use the material where there are known allergies or contact allergies.
- Do not use the material if it is not possible to prepare it in the manner described.
- Do not use the material if the prescribed machine templates for milling the blanks cannot be used.

## Basic safety information

---

- For dental use only.
- Wear mouth protection and goggles when milling, otherwise sanding dust may get into the eyes and respiratory tract.

## Notes on processing

---

- Always process the material with the specified machine templates, to avoid overheating and thus prevent damage to the material.

## Processing the blanks

---

### I. Restoration designs/design specifications

- ▶ Design the dental structures for the restoration with suitable CAD software according to the manufacturer's information.

### II. Creating the structure

**Note:** The blanks can be processed with all open CAD/CAM milling machines. Use a 5-axis milling machine to reach undercutting areas optimally.

**Note:** Milling/grinding process and machine templates must be obtained from machine manufacturers.

1. Secure the blank into the CAD/CAM milling machine.
2. Start the milling process.

### III. Removing the structure

- ▶ Remove milled structure carefully from the blank with a separator or a milling tool.

### IV. Preparing the structure

1. Carry out the rough preparation, for example with a hard metal milling tool with FSQ teeth or a black ring cutter for titanium.

**Caution!** Material damage. When polishing, keep the pressure low to prevent smearing due to heat development.

2. Carry out the fine preparation and polishing with simple, crosscut HM millers.

**Caution!** Material damage. The following are unsuitable for cleaning: Steam, autoclave, not expressly mentioned organic solvents and acids, highly concentrated cleaning alcohols, chemical prosthesis cleaners, abrasive cleaning agents.

3. Sandblast areas to be veneered (aluminum oxide, 50 – 110 µm, 2 – 3 bar).
4. Clean the structure with oil-free and water-free air.
5. Apply suitable primer (e.g. Luxatemp-Glaze & Bond) to the structure in accordance with the manufacturer's specifications. Use primers which are designed for bonding PEEK with the material used.
6. Individualize the structure as needed with an appropriate veneering resin, in accordance with the manufacturer's specifications.

## Luting restoration

---

### Temporary luting

- ▶ Lute restorations with a suitable temporary cement or glass ionomer cement in accordance with manufacturer's instructions.

## Permanent luting

**Caution!** Material damage. The following are unsuitable for cleaning: Steam, autoclave, not expressly mentioned organic solvents and acids, highly concentrated cleaning alcohols, chemical prosthesis cleaners, abrasive cleaning agents.

1. Clean the inside of the restoration with isopropyl alcohol.
2. Sandblast restoration (aluminum oxide, 50 – 110 µm, 2 – 3 bar).
3. Clean the restoration with oil-free and water-free air.
4. If necessary, clean the top of restoration with isopropyl alcohol and a clean brush.
5. Apply Luxatemp-Glaze & Bond in accordance with the manufacturer's instructions and light-cure.
6. Lute restoration with a suitable adhesive dual- or self-curing composite cement according to the manufacturer's instructions.

## Composition

Approx. 80% polyether ether ketone (PEEK), approx. 20% titanium dioxide, additives, pigments

## Storage and shelf life

- ▶ Store in a dry place at room temperature (15 – 25 °C / 59 – 77 °F).
- ▶ Protect against exposure to direct sunlight.
- ▶ Do not use after the expiration date.

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