LuxaCam

Zircon HT – Ronde

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English

Product description

LuxaCam Zircon HT – Blank is a pressed translucent milling blank from biocompatible zirconium oxide. The material is available as a blank for use in CAD/CAM technology in different colors.

Indications

- Framework structures and fully anatomical crowns and bridges in the anterior region with a maximum of 4 intermediate elements between 2 bridge abutments
- Fully anatomical crowns and bridges in the anterior region with a maximum of 2 abutments and 1 terminal bridge element
- Framework structures and fully anatomical crowns and bridges in the posterior region with a maximum of 2 intermediate elements between 2 bridge abutments
- Fully anatomical crowns and bridges in the posterior region with a maximum of 2 abutments and 1 terminal bridge element
- Veneers
- Inlays / onlays
- Primary telescopes
- Bridge constructions

Contraindications

- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use the material if it is not possible to prepare it in the manner described.
- Do not use the material if the prescribed machine templates for milling the blanks cannot be used.

Basic safety information

- For dental use only.
- Wear mouth protection and goggles when milling, otherwise sanding dust may get into the eyes and respiratory tract.

Notes on processing

- Always process the material with the specified machine templates, to avoid overheating and thus prevent damage to the material.

Processing the blanks

I. Restoration designs/design specifications

- ▶ Design the dental structures for the restoration with suitable CAD software according to the manufacturer's information. Transitions between connectors and crown/cap rounded design. For crowns and bridges, follow the general structure of the restoration.

Note: Through magnification copy milling or commercial CAD/CAM technology, blanks are milled to fully anatomical restorations. The magnification or shrink factor is specified on the blank used (see labeling on the blank).

II. Creating the structure

Note: The blanks can be processed with all open CAD/CAM milling machines. Use a 5-axis milling machine to reach undercutting areas optimally.

Note: Milling/grinding process and machine templates must be obtained from machine manufacturers.

1. Secure the blank into the CAD/CAM milling machine.
2. Start the milling process.

III. Removing the structure

- ▶ Remove milled structure carefully from the blank with a separator or a milling tool.

IV. Sintering the structure

- ▶ Sinter construction in accordance with tables.

Note: In order to obtain maximum transparency, do not use a speed program.

	Speed	Normal
Heat-up rate 1	50 °C per min to 1100 °C	10 °C per min to 950 °C
Hold time	-	
Heat-up rate 2	20 °C per min to 1500 °C	6 °C per min to 1500 °C
Hold time at final temperature	30 min	90 min
Final temperature	1500 °C	
Cooling down	Unregulated in the closed oven	

	Slow	Translucency
Heat-up rate 1	5 °C per min to 950 °C	5 °C per min to 950 °C
Hold time	-	
Heat-up rate 2	2 °C per min to 1500 °C	2 °C per min to final temperature
Hold time at final temperature	120 min	120 min
Final temperature	1500 °C	1500 °C – 1630 °C
Cooling down	Unregulated in the closed oven	

V. Preparing the structure

Caution! Material damage. Avoid overheating and thereby damage to the material in the preparation.

Caution! Material damage. Do not use sintered diamonds, corundum stones or hard metal milling. Do not sand-blast veneer surfaces.

Caution! Material damage. Do not separate with cutting discs, in order to avoid predetermined break points.

1. Adapt design as needed in the wet milling procedure with diamond milling tools (grain $\leq 30 \mu\text{m}$).
2. Carry out edge reduction and fine-tuning with fine grain diamond miller and low contact pressure.
3. Individualize design by polishing, ceramic applications and painting and glazing colors.

Note: All ceramics are suitable for zircon constructions (WAK 10 - $10.6 \cdot 10^{-6}/\text{K}$). Use ceramic in accordance with the manufacturer's instructions.

Luting restoration

Prepare the restoration

1. Clean the restoration and preparation.
2. Check occlusion and mill if necessary using suitable instruments.

Retentive luting of the restoration

Note: Pay attention to a parallel-walled preparation.

- ▶ Use a suitable cement (e.g. glass ionomer cement) in accordance with manufacturer's instructions.

Adhesive luting of the restoration

Note: Phosphoric acid can compromise the optimal bond with zircon oxide.

1. Sand-blast the parts of the restoration to be bonded (aluminum oxide, $50 \mu\text{m}$, 2 bar).
2. Clean restoration with alcohol or steam.
3. Dry the restoration with oil-free and water-free air.
4. Apply a suitable primer (e.g. LuxaBond Universal) in accordance with the manufacturer's instructions and light-cure.
5. Lute restoration with a suitable adhesive dual- or self-curing composite cement (e.g. PermaCem Universal) according to the manufacturer's instructions.

Composition

Yttrium-oxide-stabilized tetragonal zirconium dioxide

Storage and shelf life

- ▶ Store in a dry place at room temperature (15 – 25 °C / 59 – 77 °F).
- ▶ Keep away from direct sunlight.
- ▶ Do not use after the expiration date.

Packaging

1 Zircon HT blank 98 mm x 14 mm	
A1	REF 170100
A2	REF 170101
A3	REF 170102
B2	REF 170103
C2	REF 170104

1 Zircon HT blank 98 mm x 16 mm	
A1	REF 170105
A2	REF 170106
A3	REF 170107
B2	REF 170108
C2	REF 170109

1 Zircon HT blank 98 mm x 18 mm	
A1	REF 170110
A2	REF 170111
A3	REF 170112
B2	REF 170113
C2	REF 170114

1 Zircon HT blank 98 mm x 20 mm	
A1	REF 170115
A2	REF 170116
A3	REF 170117
B2	REF 170118
C2	REF 170119