# LuxaCam

## Zircon HT- Plus

DMG Chemisch-Pharmazeutische Fabrik GmbH  
Elbgaustraße 248 · 22547 Hamburg · Germany  
www.dmg-dental.com  
092015/#1\_2016-04

**DE**      Gebrauchsinformation

**EN**      Instructions for use

### Gebrauchsinformation

### Deutsch

#### Produktbeschreibung

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## English

### Product description

LuxaCam Zircon HT – Plus are pressed translucent milling blanks from biocompatible zirconium oxide. The material is available as multi-layered blank or plain block in different sizes for use in CAD/CAM technology.

### Indications

- Framework crowns and bridges with a maximum of 1 intermediate element between 2 bridge abutments
- Fully anatomical crowns and bridges with a maximum of 2 abutments and 1 terminal bridge element
- Veneers
- Inlays / onlays

### Contraindications

- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use the material if it is not possible to prepare it in the manner described.
- Do not use the material if the prescribed machine templates for milling the blanks cannot be used.

### Basic safety information

- For dental use only.
- Wear mouth protection and goggles when milling, otherwise sanding dust may get into the eyes and respiratory tract.

### Notes on processing

- Always process the material with the specified machine templates, to avoid overheating and thus prevent damage to the material.

### Processing the blanks

#### I. Restoration designs/design specifications

- ▶ Design the dental structures for the restoration with suitable CAD software according to the manufacturer's information. Transitions between connectors and crown/cap rounded design. For crowns and bridges, follow the general structure of the restoration.

**Note:** By magnification copy milling or commercial CAD/CAM technology, blanks are milled to fully anatomical restorations. The magnification or shrink factor is specified on the blank used (see labeling on the blank/block).

#### II. Creating the structure

**Note:** The blanks can be processed with all open CAD/CAM milling machines. Use a 5-axis milling machine to reach undercutting areas optimally.

**Note:** Milling/grinding process and machine templates must be obtained from machine manufacturers.

1. Secure the blank into the CAD/CAM milling machine. Follow the sequence of the colors.
2. Start the milling process.

### III. Removing the structure

- ▶ Remove milled structure carefully from the blank/block with a separator or a milling tool.

### IV. Sintering the structure

- ▶ Sinter construction in accordance with tables.

**Note:** In order to obtain maximum transparency, do not use a speed program.

	Speed	Normal
Heat-up rate 1	50 °C per min to 1100 °C	10 °C per min to 950 °C
Hold time	-	
Heat-up rate 2	20 °C per min to 1500 °C	6 °C per min to 1500 °C
Hold time at final temperature	30 min	90 min
Final temperature	1500 °C	
Cooling down	Unregulated in closed oven	

	Slow	Translucency
Heat-up rate 1	5 °C per min to 950 °C	5 °C per min to 950 °C
Hold time	-	
Heat-up rate 2	2 °C per min to 1500 °C	2 °C per min to final temperature
Hold time at final temperature	120 min	120 min
Final temperature	1500 °C	1500 °C – 1630 °C
Cooling down	Unregulated in closed oven	

### V. Preparing the structure

**Caution!** Material damage. Avoid overheating and thereby damage to the material in the preparation.

**Caution!** Material damage. Do not use sintered diamonds, corundum stones or hard metal milling. Do not sand-blast veneer surfaces.

**Caution!** Material damage. Do not separate with cutting discs, in order to avoid predetermined break points.

1. Adapt design as needed in the wet milling procedure with diamond milling tools (grain  $\leq 30 \mu\text{m}$ ).
2. Carry out edge reduction and fine-tuning with fine grain diamond miller and low contact pressure.
3. Individualize design by polishing and applying painting and glazing colors.

### Luting restoration

---

#### Prepare the restoration

1. Clean the restoration and preparation.
2. Check occlusion and grind if necessary using suitable instruments.

#### Retentive luting of the restoration

**Note:** Pay attention to a parallel-walled preparation.

- ▶ Use a suitable cement (e.g. glass ionomer cement) in accordance with manager's instructions.

#### Adhesive luting of the restoration

**Note:** Phosphoric acid can compromise the optimal bond with zircon oxide.

1. Sand-blast the parts of the restoration to be bonded (aluminum oxide, 50  $\mu\text{m}$ , 2 bar).
2. Clean restoration with alcohol or steam.
3. Dry the restoration with oil-free and water-free air.
4. Apply a suitable primer (e.g. LuxaBond Universal) in accordance with the manufacturer's instructions and light-cure.
5. Lute restoration with a suitable adhesive dual- or self-curing composite cement (e.g. PermaCem Universal) according to the manufacturer's instructions.

### Composition

---

Yttrium-oxide-stabilized tetragonal zirconium dioxide

## Storage and shelf life

- ▶ Store in a dry place at room temperature (15–25 °C / 59–77 °F).
- ▶ Keep away from direct sunlight.
- ▶ Do not use after the expiration date.

## Packaging

1 Zircon HT-Plus blank 98 mm x 14 mm		
A1		REF 170120
A2		REF 170121
A3		REF 170122
B1		REF 170123
B2		REF 170124
C2		REF 170125

1 Zircon HT-Plus blank 98 mm x 16 mm		
A1		REF 170126
A2		REF 170127
A3		REF 170128
B1		REF 170129
B2		REF 170130
C2		REF 170131

1 Zircon HT-Plus blank 98 mm x 18 mm		
A1		REF 170132
A2		REF 170133
A3		REF 170134
B1		REF 170135
B2		REF 170136
C2		REF 170137

1 Zircon HT-Plus blank 98 mm x 20 mm		
A1		REF 170138
A2		REF 170139
A3		REF 170140
B1		REF 170141
B2		REF 170142
C2		REF 170143

5 Zircon HT-Plus block 40 mm		
A1		REF 170149
A2		REF 170150
A3		REF 170151
B1		REF 170152
B2		REF 170154
C2		REF 170153