LuxaCam

Zircon LT – Ronde

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English

Product description

LuxaCam Zircon LT – Blanks are pressed translucent milling blanks from biocompatible zirconium oxide. The material is available as a blank for use in CAD/CAM technology in different sizes.

Indications

- Crown caps
- Bridge structures with up to 16 bridge elements: In the anterior region with a maximum of 4 intermediate elements between 2 bridge abutments; in the posterior region with a maximum of 2 intermediate elements between 2 bridge abutments
- Fully anatomical crowns and bridges with a maximum of 2 abutments and 1 terminal bridge element
- Primary telescopes
- Bridge constructions

Contraindications

- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use the material if it is not possible to prepare it in the manner described.
- Do not use the material if the prescribed machine templates for milling the blanks cannot be used.

Basic safety information

- For dental use only.
- Wear mouth protection and goggles when milling, otherwise sanding dust may get into the eyes and respiratory tract.

Notes on processing

- Always process the material with the specified machine templates, to avoid overheating and thus prevent damage to the material.

Processing the blanks

I. Restoration designs/design specifications

- ▶ Design the dental structures for the restoration with suitable CAD software according to the manufacturer's information. Transitions between connectors and crown/cap rounded design. For crowns and bridges, follow the general structure of the restoration.

Note: By magnification copy milling or commercial CAD/CAM technology, blanks are milled to fully anatomical restorations. The magnification or shrink factor is specified on the blank used (see labeling on the blank).

II. Creating the structure

Note: The blanks can be processed with all open CAD/CAM milling machines. Use a 5-axis milling machine to reach undercutting areas optimally.

Note: Milling/grinding process and machine templates must be obtained from machine manufacturers.

1. Secure the blank into the CAD/CAM milling machine.
2. Start the milling process.

III. Removing the structure

- ▶ Remove milled structure carefully from the blank with a separator or a milling tool.

IV. Sintering the structure

- ▶ Sinter construction in accordance with tables.

	Speed	Normal	Slow
Heat-up rate 1	50 °C per min to 1100 °C	10 °C per min to 950 °C	5 °C per min to 950 °C
Hold time	-		
Heat-up rate 2	20 °C per min to 1500 °C	6 °C per min to 1500 °C	2 °C per min to 1500 °C
Hold time at final temperature	30 min	90 min	120 min
Final temperature	1500 °C		
Cooling down	Unregulated in closed oven		

V. Preparing the structure

Caution! Material damage. Avoid overheating and thereby damage to the material in the preparation.

Caution! Material damage. Do not use sintered diamonds, corundum stones or hard metal milling. Do not sand-blast veneer surfaces.

Caution! Material damage. Do not separate with cutting discs, in order to avoid predetermined break points.

1. Adapt design as needed in the wet milling procedure with diamond milling tools (grain $\leq 100 \mu\text{m}$).
2. Carry out edge reduction and fine-tuning with fine grain diamond miller and low contact pressure.
3. Individualize design by applying ceramic.

Note: All ceramics are suitable for zircon constructions (WAK 10 - $10.6 \cdot 10^{-6}/\text{K}$). Use ceramic in accordance with the manufacturer's instructions.

Luting restoration

Prepare the restoration

1. Clean the restoration and preparation.
2. Check occlusion and grind if necessary using suitable instruments.

Retentive luting of the restoration

Note: Pay attention to a parallel-walled preparation.

- ▶ Use a suitable cement (e.g. glass ionomer cement) in accordance with manager's instructions.

Adhesive luting of the restoration

Note: Phosphoric acid can compromise the optimal bond with zircon oxide.

1. Sand-blast the parts of the restoration to be bonded (aluminum oxide, $50 \mu\text{m}$, 2 bar).
2. Clean restoration with alcohol or steam.
3. Dry the restoration with oil-free and water-free air.
4. Apply a suitable primer (e.g. LuxaBond Universal) in accordance with the manufacturer's instructions and light-cure.
5. Lute restoration with a suitable adhesive dual- or self-curing composite cement (e.g. PermaCem Universal) according to the manufacturer's instructions.

Composition

Yttrium-oxide-stabilized tetragonal zirconium dioxide

Storage and shelf life

- ▶ Store in a dry place at room temperature (15 – 25 °C / 59 – 77 °F).
- ▶ Keep away from direct sunlight.
- ▶ Do not use after the expiration date.

Packaging

1 Zircon LT blank 98 mm x 16 mm, white	REF 170144
1 Zircon LT blank 98 mm x 18 mm, white	REF 170145
1 Zircon LT blank 98 mm x 20 mm, white	REF 170146
1 Zircon LT blank 98 mm x 22 mm, white	REF 170147
1 Zircon LT blank 98 mm x 25 mm, white	REF 170148