# LuxaCore

Dual

## Inhalt | Contents

Gebrauchsinformation – DEUTSCH	4
Instructions for use – ENGLISH	14
Mode d'emploi – FRANÇAIS	23
Istruzioni d'uso – ITALIANO	33
Instrucciones de empleo – ESPAÑOL	43
Instruções de uso – PORTUGUÊS	53
Gebruiksaanwijzing – NEDERLANDS	63
Brugsanvisning – DANSK	73
Användarinstruktioner – SVENSKA	82
Instrukcja użycia – POLSKI	91
Информация по применению – РУССКИЙ	102
使用说明-简体中文	113
사용 설명서 – 독일어	121

## **Produktbeschreibung**

---

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## Product Description

---

LuxaCore-Dual is an automatic mixing, dual-cure composite that has been specially developed for all types of core build-ups and build-up fillings. Radiopaque LuxaCore-Dual stands out due to its high compressive strength, and it cuts like dentine. The curing time can be self-determined due to additional light curing. It is possible to apply LuxaCore-Dual directly using the Intraoral-Tip and the Endo-Tip.

## Indications

---

- All types of core build-up
- Cementation of root posts

## Contraindications

---

- Do not use the material in cases of allergies to any of the components or in the event of contact allergies.

- Do not use the material with a single-bottle bond because an optimum adhesive bond is generally not assured with chemical or dual-cure materials.
- Do not use the material if a dry working environment is not possible.

## Basic safety instructions

---

- ▶ Keep out of reach of children!
- ▶ For dental use only!
- ▶ In order to avoid any pulp reactions, use suitable pulp/dentine protection!
- ▶ Avoid contact with skin and eyes. In case of accidental contact, irrigate immediately with plenty of water and consult a physician, if necessary!

## Side effects

---

No systemic side effects are known.

## Interactions

---

- Materials containing eugenol may inhibit the polymerization of LuxaCore-Dual and lead to discoloration.
- Air containing water and oil can inhibit the polymerization of LuxaCore-Dual at the contact surfaces.

## Application system

---

- ▶ Automix cartridge: See “Using the Automix cartridge” on page 18.
- ▶ Smartmix syringe: See “Using the Smartmix syringe” on page 19.

## Timing

---

Working time*	≈ 1:30 min
Chemical curing	intraoral 5:00 min
Light curing	20 s (layer ≤ 2 mm) 40 s (layer ≤ 4 mm)

**\*Please note:** The times apply at a room temperature of 23 °C and a normal relative air humidity of 50 %. Higher temperatures shorten, lower temperatures increase this time.

## Instructions for use

---

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm<sup>2</sup>. Place the light as close as possible to the material.

## Recommended use

---

### Carry out core build-up

1. Isolate the tooth prior to preparation using a rubber dam.

2. Remove existing fillings and caries.
3. If the bonding agent used requires an additional etching step then etch using the etching gel per the manufacturer's instructions.

**Please note:** DMG recommends the use of the dual-cure adhesive system LuxaBond-Total Etch. In this case it is essential to etch the affected enamel and dentine surfaces. Observe the manufacturer's instructions.

**Caution!** Unprotected pulp.

- ▶ Prior to bonding in deep cavities in areas close to pulp, provide protection by using a small quantity of calcium hydroxide and a thin layer of glass ionomer cement for example.
4. Apply bonding agent in accordance with the manufacturer's instructions.
  5. In order to simplify the application of LuxaCore-Dual, place a matrix around the prepared tooth and apply LuxaCore-Dual directly into the cavity. If a pre-fabricated core build-up form is used, apply

LuxaCore-Dual directly into the core build-up form and place on the prepared tooth.

**Please note:** LuxaCore-Dual can be formed within the working time of 1:30 minutes using a standard composite instrument, e.g. a Heidemann spatula. In order to enhance processing, coat the Heidemann spatula with a bond as a separating medium.

6. Leave LuxaCore-Dual to bond intraorally for approx. 5 minutes or expose to a suitable light curing unit for a minimum 40 seconds (maximum layer thickness  $\leq 4$  mm). In the case of thinner layers ( $\leq 2$  mm) exposure to light for 20 seconds is sufficient.

**Please note:** With chemical curing, low temperatures (e.g. when working in the area of the front teeth) lead to longer bonding times.

7. If used, remove a matrix only once the material has hardened fully.
8. Produce the preparation for the desired restoration form using standard preparation instruments.

9. If a plastic temporary restoration is to be applied to the prepared tooth, isolate the preparation using a suitable separation medium (e.g. Vaseline).

### **Cementation of root posts**

1. Carry out any preparation procedures necessary for root posts in accordance with the manufacturer's instructions.
2. Prepare the root canal in accordance with the requirements of the respective root post.
3. If the bonding agent used requires an additional etching step then etch inside the canal using the etching gel per the manufacturer's instructions.

**Please note:** It is advisable to use the dual-cure adhesive system LuxaBond-Total Etch. In this case it is essential to etch the affected enamel and dentine surfaces. Observe the manufacturer's instructions.

4. Apply bonding agent in accordance with the manufacturer's instructions. If an adhesive system with a light-curing component is

used, use a paper tip to remove any excess from the canal prior to the application of light. Otherwise the lumen may be blocked.

5. Apply LuxaCore-Dual to the canal lumen. When doing so keep the end of the tip constantly immersed in the material. If desired, it is additionally possible to coat the root post with LuxaCore-Dual. The LuxaCore-Dual working time is 1:30 minutes.

**Please note:** When using LuxaBond-Total Etch, the working time of LuxaCore-Dual in the root canal is reduced to approx. 30 seconds due to the contact with the pre-bond. It is also possible that the working time may be reduced with other bonding agents. Observe the manufacturer's instructions.

6. Set the root post.
7. Using a suitable lamp, apply light to the material for 20 seconds.
8. Leave the material to harden for approx. 5:00 minutes. This guarantees complete curing,


also in the areas that cannot be reached by light.

**Please note:** With LuxaCore-Dual it is possible to work on the coronal portion of the tooth as a core build-up. Due to the easy-flowing consistency, it is recommended that work be carried out using a matrix (see “Carry out core build-up” on page 15).

## Using the Automix cartridge

---

### Insert cartridge

1. Push up the lever [A] at the back of the Automix dispenser and pull back the slide [B] completely.
2. Lift up the plastic lever [C], insert the cartridge and lock it by pushing down the plastic lever [C].

**Please note:** Make sure that the notches on the cartridge line up with those on the Automix-Dispenser.

### Attaching a mixing tip

1. Remove the cartridge cap [D] or used mixing tip by turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

**Please note:** In order to achieve an optimum mixing result, DMG recommends the use of mixing tips available from DMG. All mixing tips are also available as refill packs.

**Please note:** Make sure that the notches on the mixing tip line up with those on the cartridge.

3. Turn the mixing tip clockwise 90° to lock it into place.
4. Attach an Intraoral-Tip [F] or Endo-Tip.

### Applying the material


- ▶ When using a cartridge for the first time, express a pea-sized amount of material and discard it.
- ▶ The material is mixed in the mixing tip by pulling the trigger [G] of the Automix Dispenser and it is then applied directly.

**Please note:** After the application, leave the used Mixing-Tip on the cartridge to cap unused material!


### Remove cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] completely.
2. Lift the plastic lever [C] and remove the cartridge.

## Using the Smartmix syringe


1. Before attaching the Mixing-Tip, remove the cap [A] or used Mixing-Tip after turning it 90° counterclockwise, and discard it.


2. Attach a new Mixing-Tip [B].

**Please note:** In order to achieve an optimum mixing result, DMG recommends the use of Mixing-Tips available from DMG. All Mixing-Tips are also available as refill packs.

**Please note:** Make sure that the notches on the Mixing-Tip line up with those on the Smartmix syringe [C].

3. Turn the Mixing-Tip clockwise 90° to lock it.

4. Attach an Intraoral-Tip or Endo-Tip.
5. The material is mixed in the Mixing-Tip when it is expressed and it can be applied directly.

**Please note:** When using a syringe for the first time, express a pea-sized amount of material and discard it!

Leave the used Mixing-Tip on the Smartmix syringe after use. It serves as a cap!

## Disinfection of the Automix-Dispenser

The Automix-Dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute). The Automix Dispenser is also autoclavable.

## Technical data

Compressive strength	300 MPa
Flexural strength	100 MPa
Diametral tensile strength	50 MPa
Water absorption	20 µg/mm <sup>3</sup>

## Ingredients

Barium glass and pyrogenic silicic acid in a Bis-GMA based matrix from dental resins. Filler content: 72% by weight = 49% by vol. (0.02–2.4 µm)

## Storage and shelf life

- ▶ Store in a dry place at 2–25 °C/36–77 °F!

- ▶ If possible, store unopened material in a refrigerator.
- ▶ Do not use after the expiration date.

## Classification

ISO 4049:2000, type 2, class 3

## Available product sizes

### LuxaCore-Automix Dual

1 cartridge @ 50 g (25 ml) paste, 35 Automix-Tips, 25 Intraoral-Tips, 10 Endo-Tips*	
Shade A3	REF 212011
Shade blue	REF 212012
Shade white	REF 212126

\* In the USA and Canada, 35 Automix- & 35 Intraoral-Tips are included.

1 cartridge @ 95 g (50 ml) paste***, 25 Automix-Tips, 15 Intraoral-Tips, 10 Endo-Tips	
Shade A3	REF 212326

\*\*\* Not available in all countries.

### LuxaCore-Smartmix Dual

2 syringes @ 9 g (5 ml) paste, 20 Smartmix-Tips, 10 Intraoral-Tips, 10 Endo-Tips**	
Shade A3	REF 212039
Shade blue	REF 212038
Shade white	REF 212064

\*\* In the USA and Canada, 10 Smartmix-Tips and 10 Intraoral-Tips are included.

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**DMG** Chemisch-Pharmazeutische Fabrik GmbH