Luxatemp

Handmix

DMG Chemisch-Pharmazeutische Fabrik GmbH
Elbgaustraße 248 · 22547 Hamburg · Germany
www.dmg-dental.com
091657/#2-2015-03

CE 0482

Gebrauchsinformation

Deutsch

Produktbeschreibung

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English

Product description

Luxatemp-Handmix is a self-curing composite for the fabrication of temporary crowns and bridges, inlays, onlays and veneers. The two-component material based on multifunctional methacrylates is mixed manually.

Indications

- Fabrication of temporary crowns, bridges, inlays, onlays, partial crowns and veneers.
- Making long-term temporary prostheses.

Contra-indications

- Do not use the material in cases of allergies to any of the substances contained or, as can occur in rare cases, in the event of contact allergy.

Basic safety instructions

- Only for dental use!
- Keep out of reach of children!
- Avoid contact with skin and eyes! In the event of accidental contact, rinse thoroughly with plenty of water immediately and consult a physician if necessary.
- Wear mouth protection and goggles when constructing the temporary prosthesis, otherwise sanding dust may get into the eyes and respiratory tract!

Timing

0:00–1:00 min	Mixing and insertion into the mouth
2:30–3:00 min	Removal from the mouth
6:00–7:00 min	End of curing Work on the temporary prosthesis

Please note: The times given apply at a room temperature of 23 °C and a normal relative air humidity of 50 %. Higher temperatures shorten, lower temperatures increase these times.

Dosage and mixing

- Express the required quantities of base and catalyst paste from the tube on the mixing block in equal amounts and mix into a homogenous, bubble-free paste.

Advice on application

- Allow a quantity of paste equal to one line of the scale for each element of the bridge.
- Ensure that the thickness of the paste corresponds to the opening of the tube
- Apply the material to the impression using the enclosed application syringe.
- Close the tube again after use!

Recommended Application

1. Before preparing the crown or bridge or prior to extraction of a tooth, take an impression using alginate, silicone or a thermoplastic impression material. Alternatively, the temporary prosthesis can be made using a laboratory-prepared thermoforming matrix or, in the case of single crowns, with a strip crown. Adhere precisely to the relevant manufacturer's instructions!

Please note: in a silicone casting, undercuts will have to be evened out and if necessary drainage grooves may have to be cut in. To create a connection with the abutment teeth in the case of a gap in the molar region, cut a groove into the casting.

2. Dry the prepared teeth and even out any undercuts in the preparation. Then lightly lubricate the stumps and surrounding tissue plus any composite core build-ups using petroleum jelly or a similar separating medium.

Please note: to avoid bubbles, it is important that the syringe always remains immersed in the material during application.

3. Dispense Luxatemp-Handmix into the impression or thermoforming matrix. Then dispense the material onto the occlusal surfaces of the impression or thermoforming matrix and fill up to the gingival areas, only overbuilding slightly.
4. Place the impression or thermoforming matrix filled with Luxatemp-Handmix on the prepared teeth within at most 60 s after the start of mixing and hold in place with moderate pressure.

Please note: the setting reaction has to be monitored intra-orally (e. g. using a probe), because the curing is significantly affected by the temperature in the mouth and the temporary prosthesis can only be removed during its elastic phase.

5. The temporary prosthesis must be removed from the mouth during the elastic (rubbery) phase, i. e. during the period from 2:30 to 3:00 minutes after starting mixing.

Please note: to improve the accuracy of large provisional prostheses further, they can be put back into the impression after removal from the mouth to allow hardening to complete there.

6. Before working on it, remove the oxidation prevention layer from the surface of the temporary prosthesis (e. g. with alcohol).

Caution! When processing the temporary prosthesis, sanding dust can get into the eyes and respiratory tract.

- ▶ Wear mouth protection and goggles.
- ▶ Do not inhale the sanded dust and use an extractor.

7. From about 6:00 to 7:00 minutes after mixing, the temporary prosthesis can be processed using cross-toothed grinders, flexible disks, etc.
8. Check occlusion and grind if necessary using suitable instruments.
9. If necessary, the surface of the temporary prosthesis may be glazed with a light-curing high-gloss varnish (e. g. Luxatemp-Glaze & Bond) or polished manually with a rubber polisher or goat's hair wheel.
10. Clean the stump, dry it and thoroughly remove the separating medium.
11. Attach the temporary prosthesis.

Please note: any temporary cement can be using for fixing the temporary prosthesis. DMG recommends using automatically mixing attachment cements, e.g. TempoCemNE (non-eugenol zinc oxide cement). Cements containing eugenol can inhibit the polymerization of resin-based attachment cements during the permanent cementation.

12. Check for occlusions and adjust accordingly.

Personalization and Repair of temporary prosthesis

Please note: When personalizing and repairing a freshly-made temporary prosthesis it does not need to be roughened and a bonding agent is not required (continue with step 2).

Caution! When processing the temporary prosthesis, sanding dust can get into the eyes and respiratory tract.

- ▶ Wear mouth protection and goggles.
- ▶ Do not inhale the sanded dust and use an extractor.

1. Repairing a temporary prosthesis that is being worn: Roughen the temporary prosthesis slightly using a grinder or a sandblaster at the appropriate points. Then apply a bonding agent (e.g. Luxatemp-Glaze & Bond), in accordance with the manufacturer's instructions.

2. Apply Luxatemp-Handmix to all surfaces.
3. Press the parts together for 3:00 min.
4. After 6:00 or 7:00 minutes, the appropriate points can be finished as usual.

Please note: to speed up the polymerization, the processed temporary prosthesis can be placed in warm water at 50°C (120°F) for a few minutes. Alternatively, the temporary prosthesis can also be personalized or repaired using the light-curing, easy-flowing composite LuxaFlow Star. Follow the associated instructions for use!

Disinfecting the application syringe

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute).

Side effects

There are no known systemic side effects to date.

Composition

Glass filler materials in a matrix of multifunctional methacrylates, catalysts, stabilizers, additives. Free from methyl methacrylate and peroxides. Filler content: 43 wei.-% = 24 vol.-%. (0.02 to 1.5 µm)

Storage and shelf life

- ▶ Store at room temperature (15 to 25 °C/59 to 77 °F) in a dry place!
- ▶ Protect against exposure to direct sunlight!
- ▶ Close the tube again after use!
- ▶ Do not use after the expiration date!

Disposal of the application syringe

- ▶ Dispose of the application syringe with domestic waste after using up the pack.

Available product sizes

Starter pack

1 tube base paste @ 106 g, 1 tube catalyst paste @ 6 g, 1 mixing block, 1 spatula, 1 application syringe, colour A2	REF 110111
---	------------

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