# Luxatemp-Inlay

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### Product description

---

Luxatemp-Inlay is a light-curing composite for the temporary care of inlay preparations. Luxatemp-Inlay does not adhere to instruments when being used.

### Indications

---

- Temporary care of inlay preparations.
- Temporary obturation of implant screw access canals.

### Contra-indications

---

- ▶ Do not use the material in cases of allergies to any of the components or in the event of contact allergies.
- ▶ The placement of Luxatemp-Inlay is contraindicated when a relatively dry working area or the recommended application technique is not possible.
- ▶ Do not use directly on the exposed pulp.
- ▶ Luxatemp-Inlay is not indicated for use as temporary crown and bridge material.
- ▶ Luxatemp-Inlay should not remain in the mouth for more than 4 weeks.

### Basic safety instructions

---

- ▶ Keep out of reach of children!
- ▶ For dental use only!
- ▶ Avoid contact with skin and eyes. In case of accidental contact, irrigate immediately with plenty of water and consult a physician, if necessary!

### Side effects

---

There have been no systemic side effects to date.

### Interactions

---

- Avoid materials containing eugenol, avoid moisture and oily air; they can inhibit polymerization at the contact area.
- The use of cationic mouthwashes, plaque indicators and chlorhexidine may lead to discoloration.
- Due to similar chemical compositions, Luxatemp-Inlay may bond with resin-based relining materials.

### Timing

---

Light-cure layers with a thickness of up to 5 mm	for 20 seconds;
Light-cure layers with a thickness of up to 7 mm	for 40 seconds.

### Instructions for use

---

- ▶ When using resin-based relining materials, isolate the liner with, for example, Vaseline so that the liner is not damaged or extracted when the inlay is removed from the cavity.
- ▶ In preparations with parallel walls, Luxatemp-Inlay adheres well in the cavity. Retention can be improved by modeling Luxatemp-Inlay into the undercut of the approximal space. In this case, don't use a matrix but a wooden wedge inderdentally.
- ▶ Where minimal mechanical retention is available, the temporary may be cemented in place with a eugenol-free cement.
- ▶ When relining prefabricated crowns made of polycarbonate or methacrylate, a bonding agent such as Luxatemp-Glaze & Bond must be used for the bond between Luxatemp-Inlay and the temporary crown.
- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm<sup>2</sup>. Place the light as close as possible to the material.
- ▶ Follow the instructions for use of any other materials that are used with Luxatemp-Inlay.

### Recommended use

---

#### Temporarily fill the cavity

1. Preparation, relining and impression taking are done as usual.
2. Clean and dry the cavity.
3. Place and model a sufficient amount of Luxatemp-Inlay into the cavity using a Heidemann spatula or any other suitable instrument.
4. Light-cure layers with a thickness of up to 5 mm for 20 seconds, layers with a thickness of up to 7 mm for 40 seconds. Fill and light-cure deeper cavities in layers.
5. Remove excesses, shape and finish as usual. Polishing can be done with a polishing paste using polishing cups or discs.

#### Removal of the temporary inlay

1. Insert a suitable instrument such as a probe or a scaler into the temporary inlay and extract the material from the cavity, preferably in one piece.
2. Clean the cavity with a rubber cup/a small brush and cleaning paste.

### Storage and shelf life

---

- ▶ Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F)!
- ▶ Do not use after the expiration date.

### Composition

---

Inorganic and organic filler in a matrix of multifunctional acrylates; catalysts, stabilizers, additives.  
Free from methyl methacrylate and peroxides.

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