Silatec

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English

Product description

Pourable duplicating material based on A-Silicones.

Indications

- Fabrication of plaster or investment material duplicates

Basic safety instructions

- ▶ For dental use only!
- ▶ Keep out of reach of children!

Mixing the Material

1. Shake both Silatec bottles.
2. Squeeze Silatec out of both bottles in a 1:1 ratio and mix for at least 30 sec. until a homogenous color is achieved.

Note: Depending on the size of the model, approx. 100 g to 200 g of mixed material will be needed.

Note: When automatic mixing and dispensing units are used, shake the mixed material and then use up within 3 to 4 weeks. If applicable, fill only the amount that can be used up.

Recommended use

1. Place the model in a regular duplicating mold.

Note: The model should not be moistened and should be at room temperature (18 – 25 °C/65 – 77 °F).

2. Pour Silatec in a thin stream onto the dry model within the working time.
3. Allow Silatec to set for 30 min.
4. Remove the model by briefly blowing it out with compressed air and removing it.

Casting the duplicating mold

The mold can be cast several times using standard model cements or model plasters. When doing so, wait 30 min. before the next duplicating process. Special purpose hard plaster (ISO 6873; Type 3, e.g. Moldano, Bayer) for the cast, an ultra hard dental plaster for the forms (ISO 6873, Type 4; e.g. Fujirock II, GC) or a heat-resistant investment material (e.g. Stellavest, GC) is recommended.

Storage of the duplicating mold

- Duplicating molds made of Silatec can be stored for a long time, cast, and used several times.

Technical Data

Mixing ratio	1:1
Mixing time*	00:30 min
Processing time* (including mixing)	≤ 04:00 min
Setting time*	≥ 30:00 min
Strain in compression	≈ 16 %
Compression set	≤ 0.25 %
Reproduction accuracy	≤ 20 µm
Shore-A-hardness (1 h)	≈ 19

*The times apply at a room temperature of 23 °C and a normal relative humidity of 50%. Higher temperatures shorten, lower temperatures increase these times.

Storage and shelf-life

- ▶ Store at room temperature (15 – 25 °C/59 – 77 °F) in a dry place.
- ▶ Do not use after the expiration date.

Ingredients

Polyvinylsiloxanes, polyhydrogensiloxanes, fillers, pigments, additives, platinum catalyst

Classification

ISO 14356:2003, type 2: Duplicating Material

Available product sizes

Refill packs	
1 bottle of base @ 1 kg, 1 bottle of catalyst @ 1 kg	REF 909389

Economy pack	
1 bottle of base @ 5 kg, 1 bottle of catalyst @ 5 kg	REF 909390

Bulk pack	
1 bottle of base @ 10 kg, 1 bottle of catalyst @ 10 kg	REF 909391

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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