# Supertec

DMG Chemisch-Pharmazeutische Fabrik GmbH  
Elbgaustraße 248 · 22547 Hamburg · www.dmg-dental.com  
091760/#1\_2012-06

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Deutsch

### Produktbeschreibung

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English

### Product description

Supertec is a light-curing impression tray material used for the production of customized impression trays.

### Indication

Production of customized impression trays for precision impression taking in the field of partial and full dentures, as well as implantology

### Basic safety instructions

- ▶ For dental use only!
- ▶ Keep out of reach of children!

### Side effects

Do not apply the material in cases of allergies to any of the components or in the event of contact allergies.

### Instructions for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm<sup>2</sup>. Place the light as close as possible to the material.

### Recommended use

1. Block out undercuts on the model with a suitable material, e.g. silicone impression material.
2. Mark the tray's edge with a pencil and make a place-marker out of modeling wax.

**Note:** Exposed gypsum parts should be treated with an alginate insulation. DMG recommends that the gypsum surfaces are irrigated.

3. Remove protective cover from the prefabricated tray and adapt the tray to the model.
4. Remove excess tray material along the pencil line with a modeling knife.
5. Form a tray grip from the excess tray material and attach it to the model.
6. If necessary, smooth material junctions with petroleum jelly.
7. Cure the impression tray on the model for approx. 3 minutes with a suitable light unit, e.g. Kulzer® XS\* (UVA-VIS blended light/stroboscope).

**Note:** During curing it is possible to apply a light-curing varnish e.g. Dry Coat to the surface of the impression tray so that it is cured at the same time as the tray material. This creates a dry, shiny surface.

8. If the impression tray is not treated with light-curing varnish, remove the dispersion layer with alcohol or acetone. If light-curing varnish is used this is not necessary.
9. Let the place holder wax boil out.
10. Remove the impression tray from the model and cure the underside for approx. 3 minutes with a suitable light source e.g. Kulzer® XS\* (UVA-VIS mixed light/stroboscope).

**Note:** During curing it is possible to apply a light-curing varnish e.g. Dry Coat to the surface of the impression tray so that it is cured at the same time as the tray material. This creates a dry, shiny surface.

**Note:** The impression tray can be subsequently treated with corundum stones or fine carbide-fine cross-cut carbide burs. Wear suitable breathing protection when performing machining.

### Storage and shelf life

- ▶ Store at room temperature (15 to 25 °C/59 to 77 °F)!
- ▶ Protect against exposure to direct sunlight!
- ▶ Do not use after the expiration date.

### Composition

Glass filler and pyrogenic silicic acid in an UDMA-based matrix of dental resins.

### Technical data

Transverse strength	> 60 MPa
Transverse modulus	> 4000 MPa

### Available product sizes

50 trays upper jaw pink	REF 210410
50 trays lower jaw pink	REF 210411

\* Product names are trademarks or registered trademarks of the respective company.

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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