# DeltaFil

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English

### Product description

DeltaFil is a radio-opaque, classic glass ionomer cement for restorations. DeltaFil chemically bonds to enamel and dentin and releases fluoride ions. It is available in the shades A1, A2, A3, A3.5 and A4 based on the Vita color system.

### Intended purpose

- Permanent restorations for Class I, II: areas that are not subject to masticatory force; also suitable for the sandwich technique
- Permanent restorations for Class V; also for the sandwich technique
- Temporary restorations for Class I, II, III and V
- Extended fissure sealing
- Core build-up
- Base/liner
- Restorations in milk teeth for Class I, II, III and V

### Limitations of use

The application is contraindicated if a relatively dry working area or the use of the recommended application technique is not possible.

### Indications

- Loss of tooth substance due to caries, fractures, tooth wear or developmental disorders
- Secondary prevention of caries

### Contraindications

- Do not use directly on the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Intended users

Dentist

### Notes for use

- Clean the instruments with water immediately after use.
- Observe the manufacturer's instructions for all materials that are used with the product.
- For deep cavities, protect areas close to the pulp with a suitable liner, e.g. by using a small quantity of a firm-binding calcium hydroxide preparation.
- The time between mixing and application should be kept short to avoid premature curing of the material.
- Air bubbles should be avoided during application.
- Do not use DeltaFil without the DeltaFil Conditioner.

### Using the Handmix version

#### Powder and liquid dosing

1. Mix the glass ionomer restorative material at room temperature. The ratio of powder and liquid is 4.96 g / 1g. This corresponds to 1 level measuring spoon of powder and 1 drop of liquid.
2. Use the measuring spoon to loosen the powder in the bottle. Do not shake or turn the bottle upside down.
3. Remove the powder with the measuring spoon and level the contents of the dosing spoon on the scraper. Do not compress the powder in doing so.
4. Place a sufficient quantity of the powder onto a mixing pad.
5. Hold the bottle with the liquid vertically, press lightly and deliver the appropriate number of drops next to the powder. The drops should not contain any air bubbles.
6. Seal the powder and liquid bottles carefully after use. There should be no powder residues on the scraper. Any liquid residues on the liquid bottle should be removed.

#### Mixing

Place the powder and the liquid onto the DeltaFil mixing pad. Use a plastic spatula to divide the powder into two equal parts. Mix the first portion with the entire amount of the liquid for 15 seconds. Add the remaining powder and mix for a further 15 seconds until a homogeneous consistency has been obtained.

### Recommended use

1. Prepare the tooth as usual. The cavity should be shaped in such a way that the material layer is no thinner than 0.5 mm at any point.
2. Next, rinse the cavity thoroughly with water, removing any excess water (e.g. with a cotton pellet or by blow drying gently). The prepared surface should still be slightly damp. The dentin must not become too dry.
3. Protect the areas close to the pulp.
4. How to use DeltaFil Conditioner: Using a cotton pellet or mini applicator sponge, apply DeltaFil Conditioner to the prepared surface and allow it to set for 10 s. Then rinse thoroughly with water. Remove any excess water as described above. The surface should still look damp in preparation for the next steps.

**Note:** DeltaFil Conditioner is an aqueous polyacrylic acid solution for the treatment of the smear layer following preparation and for the conditioning of enamel and dentin for optimized bonding. It is applied as described above. Carefully reseal the DeltaFil Conditioner bottle after use.

5. Use a matrix if necessary.
6. Insert the mixed material into the cavity, pack and adapt to the margins of the cavity before modeling the surface. DeltaFil is workable for 120 s from the start of mixing.
7. During application and modeling, avoid contamination of the surface of the glass ionomer cement with saliva or water.
8. Then allow the material to cure for approx. 240 s.
9. After 360 s from the start of mixing, finish the restoration as usual, e.g. with abrasives and finishing burs. Grinding can be carried out with or without the use of water spray.

### Timing

Mixing time	30 s
Working time (23 °C, 73 °F) from the start of mixing*	120 s
Net setting time (37 °C, 98 °F)	240 s
Start time for finishing, from the start of mixing	360 s

\* The indicated times are valid for a room temperature of 23 °C/73 °F and a normal relative humidity of 50%. Higher temperatures shorten these times and lower temperatures lengthen them.

### Interactions

No interactions are currently known.

### Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

### Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material(s).
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Performance characteristics of the medical device

Complies with ISO 9917-1:2007 chemical type C and use type C

### Storage / disposal

- Store tightly sealed at 15–25 °C (59–77 °F) in a dry place protected from light!
- Do not use after the expiration date!
- Disposal must comply with national regulations!

### Composition

Ionomer glass, polyacrylic acid, water, tartaric acid, additives.

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